

**Comfort Fit-PreFab Spinal**

**Male Standard Measurements**

Size	M1	M2	M3	M4	M5	M6	M7	M8
Axilla	30" - 32"	34" - 38"	36" - 38"	38" - 40"	40" - 44"	44" - 48"	48" - 52"	52" - 56"
Waist	26" - 28"	31" - 34"	30" - 34"	34" - 36"	38" - 40"	42" - 44"	46" - 48"	49" - 52"
Trochanter	30" - 32"	35" - 39"	34" - 38"	38" - 40"	42" - 44"	46" - 48"	47" - 51"	49" - 53"
Front	14"	16.5"	17.5"	17"	18"	19"	19"	19.5"
Back	15"	17"	18.5"	18"	19"	21"	21"	21.5"

**Female Standard Measurements**

Size	F1	F2	F3	F4	F5	F6	F7	
Axilla	30" - 32"	34" - 36"	36" - 38"	39" - 41"	42" - 46"	42" - 46"	46" - 50"	
Waist	26" - 28"	30" - 32"	32" - 34"	37" - 39"	37" - 39"	42" - 44"	48" - 50"	
Trochanter	35" - 37"	34" - 36"	39" - 41"	44" - 46"	42" - 44"	47" - 49"	51" - 55"	
Front	18"	18"	18"	16"	18"	16"	17"	
Back	20"	20"	20"	18"	20"	18"	19"	

\* **Adjustable thoracic extension can be added to any model.**